

DISPOSABLE DAILY PROTECTIVE MASK

Product Description

This product has ultra-fine fiber high filtration and strong electrostatic treatment filter.


This product can filter bacteria and dust, prevent particulate dust and pollen in the air, and can also block droplets, blood, body fluids, secretions and so on.

It is light in weight, easy to store, easy to wear, breathable, comfortable, clean and hygienic.

It is widely used in factories, hospitals, electronic meters, pharmaceuticals, food, beauty salons and other places.



Three Layers Design



3rd layer:
PP Nonwoven

2nd layer:
Melt Blown Filter Pater

1st layer:
PP Nonwoven

Certification

CE / ISO 13485:2016 / EN14683:2014
Type IIR

Product Feature

Product size: 17.5 * 9.5cm 3 layers
Main ingredients: Polypropylene non-woven fabric
Filter layer: Melt-Blown filter cartridges
Nose bridge strip: Polyolefin resin
Ear hooks: nylon, spandex

Eastwest Medico Aps
Gotfred Rodes Vej 2
DK-2920 Charlottenlund, Denmark
Phone: +45 7070 1711
Mobile: +45 4027 58881-234-567-8910
Email: info@ewmedico.com

eastwest  **medico**