

PP + PE coated Isolation Gown



This isolation gown is made of 25g PP and is coated with 18g PE material to protect from liquids and fluents. It has knitted cuffs. One unit/bag and 80 units per carton.

About the manufacturer

Dymex supply disposable coveralls, lab coats, scrub suits, gowns, bouffant caps, shoe covers, face masks for chemical handling, medical usage and for biological hazards. Products are sold in Europe, North America, South America, Africa, Southeast Asia and Japan.



The factory is Quality Management System (QMS) and ISO 13485:2003 certified. Products has passed CE certification of European Union (EU) and have FDA Medical Device Listing. Dymex Surgical Masks meet EN14683 standard of EU, and also have acquired FDA 510(K).

Product Name: Impervious Isolation Gown
Product SKU: EWIG30021
Product MOQ: 100.000 normally - 500.000 at the moment
Certification: EN13795/ANSI/AAMI PB70/ISO11135/ISO11607

Unit price USD: 8.20
Lead time: 14 days after receipt of first payment
Delivery terms: DDP
Payment terms: 65% on order, 35% on delivery

Eastwest Medico Aps
 Gotfred Rodes Vej 2
 DK-2920 Charlottenlund, Denmark
 Phone: +45 7070 1711
 Mobile: +45 4027 58881-234-567-8910
 Email: info@ewmedico.com