

# DISPOSABLE FACE SHIELDS

TECHNICAL DATA SHEET



## PRODUCT INFORMATION

Disposable face shields protect healthcare workers from chemical splashes. Made from transparent PET with an anti-fog coating, they offer clear visibility. The adjustable band ensures a secure fit, covering the full face, including the sides. They can be worn with glasses or goggles and are individually packed in a plastic bag. Conforming to (EU) 2016/425, Category II standards.



## DESCRIPTION

- Anti-fog on the inside surface
- Excellent transmittance of visible light, and good stiffness, bending and punching
- PET Clear Visor: min. 0.1 mm
- Foam guard thickness 30 mm
- Panel Size: 13.19 x 9.25 inches (approx. 38.0 x 23.0 cm)
- Adjustable Strap: 40.0 x 1.5 cm (approx.)
- Latex free, fiberglass free, distortion free
- Foam headband relieves pressure visor
- OSHA-compliant

## PRODUCT DATA

<b>Type</b>	Disposable face shield
<b>Sizes</b>	One size fits all
<b>Compliance</b>	PPE Regulation (EU) 2016/425, cat II
	CE EN166
	ANSI Z87+
	AS/NZS 1337.1
<b>Materials</b>	Browguard: Foam
	Visor: PET
	Headband: Adjustable elastic strap
<b>Color</b>	Browguard: White
	Visor: Clear
	Headband: Black
<b>Visor Size</b>	19.4 x 29.8 cm
<b>Sterility</b>	Non-sterile
<b>Latex status</b>	Not made with natural rubber latex
<b>Shelf life</b>	5 years
<b>Manufacturer QMS</b>	ISO 9001



## PACKAGING METHOD: 20PCS/BAG 30BAGS/CARTON

Carton size (cm)			Total
Length	Width	Height	
44.0	48.5	54.0	600 pcs

## REORDERING

Item	Product number
Disposable Face Shield	HC SAFETY F08

## ABOUT OUR MANUFACTURER

Since 1985, Ho Cheng Safety Enterprise Co., Ltd. (HC Safety) has been a leading Taiwanese manufacturer of Personal Protective Equipment (PPE), specializing in safety glasses, earmuffs, goggles, helmets, and face shields. Their products meet international safety standards, including CE, UKCA, ANSI, and AS/NZS, with certifications from renowned laboratories such as BSI Group, DIN CERTCO, TUV, CRITT, and SGS.

HC Safety's commitment to safety, quality, and customer satisfaction has established the company as a trusted partner worldwide. They continually strive to enhance workplace safety by delivering superior PPE solutions tailored to diverse needs.

Revision 2. Date of revision: 22/05/2024.

All rights reserved for printing errors.