

DISPOSABLE PROTECTIVE APRON

TECHNICAL DATA SHEET



PRODUCT INFORMATION

This disposable sleeveless apron is designed for use in healthcare settings, offering reliable protection and comfort. Made from seamless, liquid-proof, and stain-resistant polyethylene material, it ensures durability and easy cleaning.

The apron features an integrated neckband and adjustable back tie strips for a secure and comfortable fit. Measuring 85 x 145 cm with a thickness of 50 um, it provides ample coverage for various tasks. Resistant to water and disinfectants such as ethanol (70%) and chlorine solutions (0.05% or 500 ppm), the apron is available in white and is an ideal choice for maintaining hygiene and safety in professional environments.



DESCRIPTION

Disposable PE Aprons is designed for a variety of surgical procedures where there is potential contamination from fluid penetration risk. It is made of PE material and is commonly used in emergency rooms, trauma rooms, ICU.

Subject	Value
AAMI standard	Level 3
Style	Open back
Protection	High-level fluid

PRODUCT DATA

Type	Disposable Level 3 Apron; Open Back; With protection against high-level fluids.
Sizes	Small, Medium, Large – all unisex
Compliance	AAMI Level 3 Medical Device Class 1 (MDR 2017/745)
Material	PE
Color	White
Sterility	Non-sterile
Latex status	Not made with natural rubber latex
Shelf Life	3 years
Manufacturer QMS	ISO 13485, ISO 9001



PACKAGING METHOD:

1PCS/POLYBAG 100POLYBAGS/CARTON

Carton size (cm)			Polybag (cm)			Total
Length	Width	Height	Length	Width	Height	
52.0	28.0	32.0	12.0	19.0	3.0	100 pcs



REORDERING

Item	Product number
Disposable plastic aprons	IL0023LBTU-UN

ABOUT OUR MANUFACTURER

Founded in 1989, Medtecs International Corporation Limited is a listed company in Taiwan and Singapore, and a premier supplier and distributor of personal protective equipment (PPE) and a provider of logistics services to healthcare institutions. With a strong presence across the United States, Europe, and the Asia Pacific, Medtecs operates 14 manufacturing facilities in the Philippines, Cambodia, and Taiwan, ensuring stable production capacity.

Medtecs offers a comprehensive range of products, including medical uniforms, apparel, and consumer health items, all designed to meet diverse needs and daily usages. Their products conform to international standards, such as CE and FDA, and undergo rigorous testing to ensure safety and optimal effectiveness.

With over three decades of experience, Medtecs has built an extensive track record serving major healthcare institutions in Taiwan, the United States, Singapore, and the Philippines. Their commitment to quality, innovation, and responsible manufacturing has established them as a trusted partner in the global healthcare industry.

Revision 5. Date of revision: 23/06/2024.

All rights reserved for printing errors.